

Flow Out Viscometer

For quick testing of the viscosity of sour milk products, such as yogurt, sour milk, sour cream and cream . Simply fill the funnel with 100 ml, let the sample drain through the outlet valve and measure the time.

Included in the delivery:

- 08.3001 Flow out viscometer with nozzle Ø 2 and 3 mm
- 08.30003 Stand garnish
- Pharmacist patel 130 mm
- Pharmacist patel 200 mm
- Ring for flow out viscometer

Order No.	Type	Qty per unit
08.3000	Flow out viscometer complete	1



Accessories for Flow Out Viscometer

Order No.	Type	Qty per unit
08.3002	Nozzle for viscometer Ø 2 mm	1
08.3003	Nozzle for viscometer Ø 3 mm	1
08.3004	Nozzle for viscometer Ø 4 mm	1
08.3005	Nozzle for viscometer Ø 5 mm	1
08.3006	Nozzle for viscometer Ø 6 mm	1
08.3008	Nozzle for viscometer Ø 8 mm	1
9.262410	Stopwatch Addition stopper Amigo, impact-resistant, unbreakable main spring, ABS case, Ø 55 mm, with lanyard (120 cm), division sec: 1/10, display min: 15	1
9.260189	Stopwatch Digi clock timer TR 119 OS, large 6-digit LCD display, programmable for up to 23 h, 59 min, 59 sec 1-minute audible alarm when timed period has expired, countdown, countup, time display, clip and stand, instructions in 3 languages, supplied with 1 battery AAA 1.5 VAC	1

Glass Viscometer

Order No.	Type	Qty per unit
08.3009	100 ml with integrated nozzle Ø 2.8 mm, made of DURAN glass	1