

## Special examinations of milk- and other powders

### Determination of the Jolted Specific Weight

#### Jolting Volumeter STAV II

To determine of jolted volume according to ISO, ASTM and DIN standards.

Included in the delivery:

- Jolting volumeter STAV II
- Graduated measuring cylinder 250 ml, made of glass

Order No.	Type	Qty per unit
09.0200	Jolting volumeter STAV II complete, 230 VAC/50 Hz	1
09.0200-115	Jolting volumeter STAV II complete, 115 VAC/60 Hz	1



#### Accessories for STAV II

Order No.	Type	Qty per unit
09.0201	Spare glass cylinder 250 ml, length 330 mm	1
09.0202	Acrylic glass cover incl. frame for noise damping hood	1
09.0203	Preselection counter	1
09.0210	Noise damping hood KU H550 for measuring cylinder 10 to 250 ml	1

#### Jolting Volumeter Gerber

To determine of specific weight of a jolted sample.  
Automatic apparatus with up to 100 strokes in 30 sec.

Included in the delivery:

- 09.0220 Jolting volumeter
- 09.02201 Cylinder, 500 ml, made of stainless steel
- Country-specific power cord

Order No.	Type	Qty per unit
09.0220	Jolting volumeter Gerber complete, 230 VAC/50 Hz	1
09.0221	Jolting volumeter Gerber complete, 115 VAC/60 Hz	1

