

## Special examinations of milk- and other powders

### Determination of Miscibility and Clumping



#### Mixer for the Determination of Miscibility and Clumping

To mix powders and liquids.

3 fixed programs:

- Program 1: Speed 50 rpm, time 180 sec
- Program 2: Speed 60 rpm, time 30 sec
- Program 3: Speed 60 rpm, time 15 sec

1 variable program:

- Speed 40 to 80 rpm, time 15 to 300 sec

Specifications:

- Digital display
- Indicator when changing from fixed to variable program
- Low-noise revolve
- Power: 100 W
- Voltage: 230 VAC/50 Hz

Order No.	Type	Qty per unit
09.0300	Mixer with drum and holder without sample tube	1




09.0301+09.0303

#### Accessories for Mixer for the Determination of Miscibility

Order No.	Type	Qty per unit
09.0301	Tripod for sample tube (09.0303), made of PVC	1
09.0302	Funnel, made of stainless steel, Ø 120/36 x 185/110 mm	1
09.0303	Sample tube, Ø 52 x 245 mm, RN 45/40	1
9.230145	Stopper, hexagonal, RN 45/40, PE, without drillings	1

### Determination of Flowability of Powders



#### Flowability Funnels

Order No.	Type	Qty per unit
09.0400	Funnel battery made of stainless steel dimensions (WxDxH): 80 x 20 x 20 cm, with 7 funnels Ø 9 cm apertures Ø : 2.5, 3.0, 3.5, 4.0, 4.5 and 5.0 mm	1