

Filtration | stainless steel

Filter Ø47mm, 100ml / 500ml stainless steel filter holder

Ideal for:

- filtration requiring the sterilization such as microbiology test or use for purification of corrosive liquid, solvent, etc.




VF2 - silicone stopper
Ø47mm, 100ml



VF11 - silicone stopper
Ø47mm, 500ml


SPIN-LOCK CONNECTIONS



VF2

| VF2 - 100ml funnel | | | | | |
|--------------------|--|--------------------|---------|-----------------|-----------------|
| # | Item | Material | Size | Order number | |
| | 47mm Stainless steel filtration set, complete | | | | 00829 167200-02 |
| | 47mm Stainless steel filtration set, filter holder | | | | 00829 167103-10 |
| | Filter holder | Stainless steel | Ø47mm | | |
| 1 | Funnel | Stainless steel | 100 ml | 00829 167103-61 | |
| 2 | Support screen | Stainless steel | | 00829 167100-31 | |
| 3 | Support base | Stainless steel | | 00829 167103-60 | |
| 4 | Stopper | Silicone | No. 8 | 00829 167110-16 | |
| 5 | Receiver flask | Borosilicate glass | 1000 ml | 00829 167300-08 | |
| 6 | Fixing sucker | Silicone | | 00829 167200-45 | |

| VF11 - 500ml funnel | | | | | |
|---------------------|--|--------------------|---------|-----------------|-----------------|
| # | Item | Material | Size | Order number | |
| | 47mm Stainless steel filtration set, complete | | | | 00829 167200-11 |
| | 47mm Stainless steel filtration set, filter holder | | | | 00829 180100-50 |
| | Filter holder | Stainless steel | Ø47mm | | |
| 1 | Funnel | Stainless steel | 500 ml | 00829 167103-65 | |
| 2 | Support screen | Stainless steel | | 00829 167100-31 | |
| 3 | Support base | Stainless steel | | 00829 167103-60 | |
| 4 | Stopper | Silicone | No. 18 | 00829 167110-18 | |
| 5 | Receiver flask | Borosilicate glass | 2000 ml | 00829 167200-13 | |



VF11

Filtration | stainless steel

Filter Ø47mm, 100ml / 300ml stainless steel filter holder

Features:

- spin-lock connection facilitates fast and stable installation without clamp
- stainless steel funnel offers excellent chemical resistance and can be sterilized
- built-in overflow protection to keep the filtrate from overflowing when waste bottle is full




LF 32 - PC receiver flask
Ø47mm, 100ml



LF 33 - PC receiver flask
Ø47mm, 300ml


SPIN-LOCK CONNECTIONS



LF 32

| LF 32 - 100ml funnel | | | | | | |
|----------------------|--|-----------------|---------|--------------------|--------------------|--|
| # | Item | Material | Size | Order number | | |
| | 47mm stainless steel filtration set, complete | | | | 00829 197010-32 | |
| | 47mm stainless steel filtration set, filter holder | | | | 00829 167103-10 | |
| 1 | Filter | Stainless steel | Ø47mm | | | |
| 2 | Funnel | Stainless steel | 100 ml | 00829 167103-61 | | |
| 3 | Support screen | Stainless steel | | 00829 167100-31 | | |
| 4 | Support base | Stainless steel | | 00829 167103-60 | | |
| 5 | Stopper | Silicone | No. 16 | 00829 167110-17 | | |
| | PC Waste bottle with cover | | | | 00829 197000-11-PC | |
| 6 | Cap | PP | | 00829 197000-40 | | |
| 7 | Cover | PC | | 00829 197000-32-PC | | |
| | Overflow protection | | | | 00829 197000-39 | |
| 8 | Waste bottle | PC | 1200 ml | 00829 197000-31-PC | | |

| LF 33 - 300ml funnel | | | | | | |
|----------------------|--|-----------------|---------|--------------------|--------------------|--|
| # | Item | Material | Size | Order number | | |
| | 47mm stainless steel filtration set, complete | | | | 00829 197010-33 | |
| | 47mm stainless steel filtration set, filter holder | | | | 00829 167103-30 | |
| 1 | Filter | Stainless steel | Ø47mm | | | |
| 2 | Funnel | Stainless steel | 300 ml | 00829 167103-63 | | |
| 3 | Support screen | Stainless steel | | 00829 167100-31 | | |
| 4 | Support base | Stainless steel | | 00829 167103-60 | | |
| 5 | Stopper | Silicone | No. 16 | 00829 167110-17 | | |
| | PC Waste bottle with cover | | | | 00829 197000-11-PC | |
| 6 | Cap | PP | | 00829 197000-40 | | |
| 7 | Cover | PC | | 00829 197000-32-PC | | |
| | Overflow protection | | | | 00829 197000-39 | |
| 8 | Waste bottle | PC | 1200 ml | 00829 197000-31-PC | | |



LF 33

Filter holder | stainless steel


Filter Ø47mm, 100ml - 500ml stainless steel filter holder

Features:

- spin-lock connection facilitates fast and stable installation without clamp
- stainless steel funnel offers excellent chemical resistance and can be sterilized



SPIN-LOCK CONNECTIONS



| 100ml / 300ml / 500ml stainless steel filter holder | | | | |
|---|---|-----------------|--------------------------------------|---|
| # | Item | Material | Size | Order number |
| | 100ml, 47mm stainless steel filter holder, complete | | | 00829 167103-10 |
| | 300ml, 47mm stainless steel filter holder, complete | | | 00829 167103-30 |
| | 500ml, 47mm stainless steel filter holder, complete | | | 00829 180100-50 |
| 1 | Funnel | Stainless steel | 100 ml 300 ml 500 ml | 00829 167103-61 00829 167103-63 00829 167103-65 |
| 2 | Support screen | Stainless steel | | 00829 167100-31 |
| 3 | Support base | Stainless steel | | 00829 167103-60 |
| 4 | Stopper | Silicone | No. 16 (100/300ml) No. 18 (500ml) | 00829 167110-17 00829 167110-18 |

| Stainless steel lids for filter holder | |
|---|-----------------|
| Item | Order number |
| Stainless steel lid for 100ml filter holder | 00829 167103-11 |
| Stainless steel lid for 300ml filter holder | 00829 167103-13 |
| Stainless steel lid for 500ml filter holder | 00829 167103-15 |