

Special examinations of milk- and other powders

Determination of Wettability of Powders

Funnel Battery

For standard measurement of the time, that the powder needs, to absorb deionized water.

Order No.	Type	Qty per unit
09.0450	Funnel battery made of stainless steel dimensions (WxDxH): 50 x 10 x 15.5/20 cm, with 4 funnels Ø 5 cm	1
09.0451	Powder dosing cup with scraper	1



Determining the Angle of Repose of Powders and Granulates

Pourability Tester

Order No.	Type	Qty per unit
09.0460	Standard version according to DIN/ISO 4324 Included in the delivery: - 09.0462 Pourability Tester - 09.04600 Funnel with outlet Ø 10 mm - 09.04601 Acrylic glass pane Ø 100 x 25 mm	1
09.0461	Modified version for milk powders with high fat content Included in the delivery: - 09.0462 Pourability Tester - 09.04610 Funnel with outlet Ø 12 mm - 09.04611 Acrylic glass pane Ø 100 x 15 mm	1



Determination of Color

V-Samples Scraper

Order No.	Type	Qty per unit
09.8135	V-samples scraper, small model Included in the delivery: - 09.81350 V-samples scraper, small - 09.8136 Sample cup black Ø 90 mm - Laboratory dish (WxDxH) 310 x 370 x 75 mm	1
09.8138	V-samples scraper, large model Included in the delivery: - 09.81380 V-samples scraper, large - 09.8139 Sample cup black Ø 142 mm - Laboratory dish (WxDxH) 310 x 370 x 75 mm	1

